

Product datasheet for **RC202568L4V**

RBM19 (NM_016196) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RBM19 (NM_016196) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RBM19
Synonyms:	Mrd1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016196
ORF Size:	2880 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202568).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016196.1
RefSeq Size:	4217 bp
RefSeq ORF:	2883 bp
Locus ID:	9904
UniProt ID:	Q9Y4C8
Cytogenetics:	12q24.13-q24.21
Domains:	RRM, RRM_1
MW:	107.3 kDa


[View online »](#)

Gene Summary:

This gene encodes a nucleolar protein that contains six RNA-binding motifs. The encoded protein may be involved in regulating ribosome biogenesis. Multiple alternatively spliced variants, encoding the same protein, have been identified.[provided by RefSeq, Apr 2009]