

## Product datasheet for **RC202574**

### NDFIP1 (NM\_030571) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	NDFIP1 (NM_030571) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NDFIP1
Synonyms:	N4WBP5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC202574 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGTTGGCGTTGGCGGCGCTGGCGGCGGTCGAGCCGGCTGCGGCAGCCGGTACCAGCAGTTGCAGA  
ATGAAGAAGAGTCTGGAGAACCTGAACAGGCTGCAGGTGATGCTCCTCCACCTTACAGCAGCATTCTGC  
AGAGAGCGCAGCATATTTTGACTACAAGGATGAGTCTGGGTTTCAAAGCCCCATCTTACAATGTAGCT  
ACAACACTGCCAGTTATGATGAAGCGGAGAGGACCAAGGCTGAAGCTACTATCCCTTTGGTTCTGGGA  
GAGATGAGGATTTGTGGGTCGGGATGATTTTGATGATGCTGACCAGCTGAGGATAGGAAATGATGGGAT  
TTTCATGTAACTTTTTTCATGGCATTCTCTTAACTGGATTGGGTTTTTCTGTCTTTTTGCCTGACC  
ACTTCAGCTGCAGGAAGGTATGGGGCCATTTCAAGATTTGGTCTCTCTAATTAATGGATCCTGATTG  
TCAGGTTTTCCACCTATTTCCCTGGATATTTGATGGTCAGTACTGGCTCTGGTGGGTGTTCTTTGTTTT  
AGGCTTCTCCTGTTTCTCAGAGGATTTATCAATTATGCAAAAGTTCGGAAGATGCCAGAACTTTCTCA  
AATCTCCCAGGACCAGAGTTCTCTTTATTTAT

**ACGCGT**ACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC202574 protein sequence  
Red=Cloning site Green=Tags(s)


MALALAALAAVEPACGSRYQQLQNEEESGEPEQAAGDAPPPYSSISAESAAYFDYKDESGFPKPPSYNVA  
 TTLPSYDEAERTKAEAT IPLVGRDEDFVGRDDFDDADQLRIGNDGIFMLTFFMAFLFNWIGFLLSFCLT  
 TSAAGRYGAISGFGLSLIKWILIVRFSTYFPGYFDGQYWLWWVFLVLGFLFLRGRFINYAKVRKMPETFS  
 NLPRTVLFYIY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk6233\\_g03.zip](https://cdn.origene.com/chromatograms/mk6233_g03.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_030571

**ORF Size:** 663 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_030571.4](#)

**RefSeq Size:** 3599 bp

**RefSeq ORF:** 666 bp

**Locus ID:** 80762

**UniProt ID:** [Q9BT67](#)


**Cytogenetics:** 5q31.3

**Protein Families:** Transmembrane


**MW:** 24.9 kDa

**Gene Summary:** The protein encoded by this gene belongs to a small group of evolutionarily conserved proteins with three transmembrane domains. It is a potential target for ubiquitination by the Nedd4 family of proteins. This protein is thought to be part of a family of integral Golgi membrane proteins. [provided by RefSeq, Jul 2008]

### Product images:



Circular map for RC202574



Western blot validation of overexpression lysate (Cat# [LY410830]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC202574 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).