

## Product datasheet for **RC202577L2V**

### SPX (NM\_030572) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SPX (NM_030572) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPX
Synonyms:	C12orf39
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_030572
ORF Size:	348 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202577).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_030572.1</a>
RefSeq Size:	638 bp
RefSeq ORF:	351 bp
Locus ID:	80763
UniProt ID:	<a href="#">Q9BT56</a>
Cytogenetics:	12p12.1
Protein Families:	Secreted Protein, Transmembrane
MW:	13.3 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a hormone involved in modulation of cardiovascular and renal function. It has also been shown in rats to cause weight loss. Several transcript variants have been found for this gene. [provided by RefSeq, Feb 2016]