

## Product datasheet for **RC202585L3V**

### Apolipoprotein L 2 (APOL2) (NM\_030882) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Apolipoprotein L 2 (APOL2) (NM_030882) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Apolipoprotein L 2
Synonyms:	APOL-II; APOL3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_030882
ORF Size:	1011 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202585).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_030882.2</a>
RefSeq Size:	2545 bp
RefSeq ORF:	1014 bp
Locus ID:	23780
UniProt ID:	<a href="#">Q9BQE5</a>
Cytogenetics:	22q12.3
Protein Families:	Druggable Genome, Transmembrane
MW:	37.1 kDa



[View online »](#)

**Gene Summary:**

This gene is a member of the apolipoprotein L gene family. The encoded protein is found in the cytoplasm, where it may affect the movement of lipids or allow the binding of lipids to organelles. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]