

Product datasheet for RC202593L2V

SH3GLB2 (NM_020145) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SH3GLB2 (NM_020145) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SH3GLB2
Synonyms:	PP6569; PP9455; RRIG1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_020145
ORF Size:	1185 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202593).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020145.2 , NP_064530.1
RefSeq Size:	2296 bp
RefSeq ORF:	1188 bp
Locus ID:	56904
UniProt ID:	Q9NR46
Cytogenetics:	9q34.11
Domains:	SH3, BAR
Protein Pathways:	Endocytosis


[View online »](#)

MW: 44 kDa