

Product datasheet for **RC202599A1V**

Human TRADD (NM_003789) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TRADD (NM_003789) AAV Particle
Tag:	Myc-DDK
Symbol:	TRADD
Synonyms:	Hs.89862
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202599 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAGCTGGGCAAAATGGGCACGAAGAGTGGGTGGGCAGCGCATACCTGTTTGTGGAGTCCTCGCTGG
ACAAGGTGGTCCTGTTCGATGCCTACGCGCACCCAGCAGAAGGTGGCAGTGTACAGGGCTCTGCAGGC
TGCTTTGGCAGAGAGCGCGGGAGCCCGACGTGCTGCAGATGCTGAAGATCCACCGCAGCGACCCGAG
CTGATCGTGCAGCTGCGATTCTCGGGCGGCAGCCCTGTGGCCGCTTCCTCCGCGCTACCGCGAGGGGG
CGCTGCGCGCCGCGCTGCAGAGGAGCCTGGCGCCGCGCTCGCCAGCACTCGGTGCCGCTGCAACTGGA
GCTGCGCGCCGCGCGAGCGGCTGGACGCTTTGCTGGCGGACGAGGAGCGCTGTTGAGTTGCATCCTA
GCCCAGCAGCCCGACCGCTCCGGGATGAAGAACTGGCTGAGCTGGAGGATGCGCTGCGAAATCTGAAGT
GCGGCTCGGGGGCCCGGGTGGCGACGGGGAGGTGCTTCGGCCCCCTTCAGCCCCCGGTGCCCTCTCT
GTCGGAGGTGAAGCCGCGCCGCGCCACCTGCCAGACTTTTCTGTTCCAGGGTCAGCCTGTAGTG
AATCGGCCGCTGAGCCTGAAGGACCAACAGACGTTGCGCGCTCTGTGGTCTCAAATGGCGCAAGGTGG
GGCGTCACTGCAGCGAGGCTGCCGGGCGCTGCGGGACCCGGCGCTGGACTCGTGGCTACGAGTACGA
GCGCGAGGAGTGTACGAGCAGGCTTCAGCTGCTGCGGCGCTTCGTGAGGCCGAGGGCCGCCGCC
ACGCTGCAGCGCTGGTGGAGGCACTCGAGGAGAAGAGCTCACCAGCCTGGCAGAGGACTTGCTGGGCC
TGACCGATCCCAATGGCGGCCTGGCC

ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202599 protein sequence Red=Cloning site Green=Tags(s) MAAGQNGHEEWVGSAYLFVSSLDKVVLSDAYAHPQQKVAVYRALQAALAESGGSPDVLQMLKIHRSDPQ LIVQLRFCGRQPCGRFLRAYREGALRAALQRSLAAALAQHSVPLQLELRAGAERLDALLADEERCLSCIL AQQPDRLRDEELAELEDALRNLCGSGARGGDGEVASAPLQPPVPSLSEVKPPPPPPPAQTFLFQGQPVV NRPLSLKDQQTFAFSVGLKWRKVGRSLQRCRALRDPALDSLAYEYEREGLYEQAFQLLRRFVQAEGRRA TLQRLVEALEENELTSLAEDLLGLTDPNGGLA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_003789
ORF Size:	936 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003789.3</u>
RefSeq Size:	1496 bp
RefSeq ORF:	939 bp
Locus ID:	8717
UniProt ID:	<u>Q15628</u>
Cytogenetics:	16q22.1
MW:	34.2 kDa