

Product datasheet for **RC202601A1V**

Human METTL7A (NM_014033) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human METTL7A (NM_014033) AAV Particle
Tag:	Myc-DDK
Symbol:	METTL7A
Synonyms:	AAM-B; AAMB
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202601 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGCTTACCATCTTTATCCTGAGACTGGCCATTTACATCCTGACATTTCCCTTGACCTGCTGAACT
TTCTGGGCTTGTTGGAGCTGGATATGCAAAAAATGGTCCCTACTTCTTGGTGAGGTTCACTGTGATATA
CAACGAACAGATGGCAAGCAAGAAGCGGGAGCTCTTCAGTAACCTGCAGGAGTTTGGGGCCCCCTCCGGG
AAACTCTCCCTGCTGGAAGTGGGCTGTGGCACGGGGGCCAACTTCAAGTTCTACCCACCTGGGTGCAGGG
TGACCTGTATTGACCCCAACCCCACTTTGAGAAGTTTTGATCAAGAGCATTGCAGAGAACCGACACCT
GCAGTTTGAGCGCTTTGTGGTAGCTGCCGGGGAGAACATGCACCAGGTGGCTGATGGCTCTGTGGATGTG
GTGGTCTGCACCCTGGTGTGCTCTGTGAAGAACCAGGAGCGGATTCTCCGCGAGGTGTGCAGAGTGC
TGAGACCGGGAGGGGCTTTCTATTTTCATGGAGCATGTGGCAGCTGAGTGTTCGACTTGGAATTACTTCTG
GCAACAAGTCCTGGATCCTGCCTGGCACCTTCTGTTTGATGGGTGCAACCTGACCAGAGAGAGCTGGAAG
GCCCTGGAGCGGGCCAGCTTCTTAAGCTGAAGCTGCAGCACATCCAGGCCCACTGTCTGGGAGTTGG
TGCGCCCTCATATCTATGGATATGCTGTGAAA

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202601 protein sequence Red=Cloning site Green=Tags(s) MELTIFILRLAIYILTFPLYLLNFLGLWSWICKKWFYFLVRFTVIYNEQMASKKRELF SNLQEFAGPSG KLSLLEVGCGTGANFKFYPPGCRVTCIDPNPNFEKFLIKSIAENRHLQFERFVVAAGENMHQVADGSVDV VVCTLVLCVKNQERILREVCRLRPGGAFYFMEHVAAECSTWNYFWQQVLDPAWHLDFDGCNL TRESWK ALERASF SKLKLQHIQAPLSWELVRPHIYGAVK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_014033
ORF Size:	732 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014033.3</u>
RefSeq Size:	3390 bp
RefSeq ORF:	735 bp
Locus ID:	25840
UniProt ID:	<u>Q9H8H3</u>
Cytogenetics:	12q13.12
MW:	28.3 kDa