

Product datasheet for **RC202615A1V**

Human FBXO27 (NM_178820) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FBXO27 (NM_178820) AAV Particle
Tag:	Myc-DDK
Symbol:	FBXO27
Synonyms:	FBG5; Fbx27
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202615 representing NM_178820 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGGCGCCTCGGTCTCCAGGGGCCGGGCCGCCGGGTCCCCGCGCCGGAGCCGGAACCCGAAGAGGCGC
TGGACCTGAGCCAACACCCCGAGAGCTGCTTCTGGTGGTGCTGAGCCACGTCCCCCGCGCACGCTGCT
CGGGCGCTGCCCAAGTGTGCCGGGGCTGGCGAGCCCTGGTGGACGGCCAGGCCCTGTGGCTGCTGATC
CTGGCCCGCGACACGGCGCCACCGGCCGCGCGCTGCTGCACCTCGCCCGCAGCTGCCAGTCTCCCGCCC
GTAACGCCAGGCCTTGCCCCCTGGGCCGCTTCTGCGCGCGCAGACCCATCGGACGCAACCTTATTCGCAA
CCCCTGCGGCCAAGAAGGCCTCCGAAAGTGGATGGTGCAACACGGTGGGACGGCTGGGTGGTGGAGGAA
AACAGGACAACCGTGCCTGGGGCCCCCTCTCAGACGTGCTTCGTGACTTCATTACAGCTGGTGTGCAAGA
AGCAGGTCTTGACCTAGAGGAGGAGGGTCTGTGGCCAGAAGTCTGGATAGTGGCAGGATTGAGATTTG
TGTCTCTGACTGGTGGGAGCCCGACACGACAGCGGCTGTATGTACAGACTCCTCGTCCAACCTTAGAC
GCCAACCAGACTGTTCTAGATAAATTCTCTGCTGTGCCTGATCCCATCCCGCAGTGGAACAACAATGCCT
GCCTTCACGTACCCACGTGTTCTCAACATCAAGATGGGCGTCCGCTTTGTGTCTTTCGAACACCGGGG
CCAGGACACACAGTTCTGGGCTGGCCACTATGGAGCCCGTGTGACCAACTCCAGTGTGATCGTGGCAGTC
CGTCTGTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202615 representing NM_178820 Red=Cloning site Green=Tags(s) MGASVSRGRAARVPAPEPEPEEALDLSQLPELLLVL SHVPPRTLGRGRQVCRGWALVDGQALWLLI LARDHGATGRALLHLARSCQSPARNARPCPLGRFCARRPIGRNLIRNPGQEGLRKWMVQHGGDGWVVEE NRTTVPGAPSQTCFVTSFSWCCKKQVLDLEEEGLWPELLDSGRIEICVSDWWGARHDSGCMYRLLVQLLD ANQTVLDKFSAVPDPIQWNNACLHVTHVFSNIKMGVRFVSFEHRGQDTQFWAGHYGARVTNSSVIVRV RLS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_178820
ORF Size:	849 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_178820.3</u>
RefSeq Size:	2338 bp
RefSeq ORF:	852 bp
Locus ID:	126433
UniProt ID:	<u>Q8NI29</u>
Cytogenetics:	19q13.2
MW:	31.4 kDa