

## Product datasheet for **RC202619L1V**

### **GJA5 (NM\_005266) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	GJA5 (NM_005266) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GJA5
Synonyms:	ATFB11; CX40
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005266
ORF Size:	1074 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202619).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005266.5</a>
RefSeq Size:	3201 bp
RefSeq ORF:	1077 bp
Locus ID:	2702
UniProt ID:	<a href="#">P36382</a>
Cytogenetics:	1q21.2
Domains:	CNX
Protein Families:	Transmembrane



[View online »](#)

**MW:** 40.4 kDa

**Gene Summary:** This gene is a member of the connexin gene family. The encoded protein is a component of gap junctions, which are composed of arrays of intercellular channels that provide a route for the diffusion of low molecular weight materials from cell to cell. Mutations in this gene may be associated with atrial fibrillation. Alternatively spliced transcript variants encoding the same isoform have been described. [provided by RefSeq, Jul 2008]