

Product datasheet for **RC202629A1V**

Human PRR15L (NM_024320) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PRR15L (NM_024320) AAV Particle
Tag: Myc-DDK
Symbol: PRR15L
Synonyms: ATAD4
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202629 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACGACTGAAATTGGTTGGTGAAGCTGACTTTCCTCCGAAAAAGAAATCCACTCCCAAAGTGCTGT
ATGAGATCCCTGACACCTATGCCAAACAGAGGGAGATGCAGAACCCCGAGGCCTGACGCTGGAGGCC
CAACAGCGACTTTAACACCCGCCTGGAGAAGATTGTGGACAAGAGCACAAAGGCAAGCACGTCAGGTCTC
TCCAACCTCAGGACGCTTCAAGGAGAAGAAGAAAGTGAGAGCCACGCTGGCAGAGAACCCTAACCTCTTTG
ATGATCACGAGGAAGGACGGTCATCAAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202629 protein sequence
Red=Cloning site Green=Tags(s)

MTTEIGWWKLTFLRKKKSTPKVLYEIPDTYAQTEGDAEPPRPDAGGPNSDFNTRLEKIVDKSTKGKHVKV
SNSGRFKEKKKVRATLAENPNLFDDHEEGRSSK

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_024320
ORF Size: 309 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_024320.2
RefSeq Size:	1530 bp
RefSeq ORF:	312 bp
Locus ID:	79170
UniProt ID:	Q9BU68
Cytogenetics:	17q21.32
MW:	11.7 kDa