

Product datasheet for **RC202633L4V**

ABHD3 (NM_138340) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ABHD3 (NM_138340) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABHD3
Synonyms:	LABH3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138340
ORF Size:	1227 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202633).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138340.3
RefSeq Size:	2085 bp
RefSeq ORF:	1230 bp
Locus ID:	171586
UniProt ID:	Q8WU67
Cytogenetics:	18q11.2
Protein Families:	Transmembrane
MW:	46 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein containing an alpha/beta hydrolase fold, which is a catalytic domain found in a very wide range of enzymes. The function of this protein has not been determined. [provided by RefSeq, Jul 2008]