

Product datasheet for **RC202645A1V**

Human BRMS1L (NM_032352) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human BRMS1L (NM_032352) AAV Particle
Tag:	Myc-DDK
Symbol:	BRMS1L
Synonyms:	BRMS1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202645 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCCAGTCCATTCCCGAGGGGATAAGAAGGAGACCAACCATCACGATGAGATGGAGGTGGACTACGCCG
AAAAATGAGGGGAGCAGCTCCGAGGACGAGGACACTGAGAGCTCGTCGGTCTCCGAGGATGGAGATAGCTC
AGAAATGGATGATGAAGACTGTGAAAGAAGAAGAATGGAATGTTTGGATGAAATGTCCAATCTTGAAAAA
CAGTTTACCGATCTCAAAGATCACTTTATAAAGAACGATTAAGTCAGGTGGATGCAAACTACAAGAAG
TCATAGCTGAAAAAGCACCAGAATACTTGGAACCGCTGGCAACTTTACAGGAAAAATGCAAAATTCGTAC
AAAGGTAGCAGGAATCTATAGAGAGCTCTGCTTAGAATCTGTAAAGAACAATATGAATGTGAAATTCAA
GCTTCTCGCCAGCATTGTGAGAGCGAAAAGCTGTTGCTATATGATACAGTCCAGAGTGAAGTACAGGAGA
AGATAAGAAGGCTTGAAGAGGATAGGCACAGCATTGATATTACCTCAGAGCTGTGGAATGATGAGCTTCA
GTCAAGAAAAAGAGGAAGGATCCTTTTCAGTCTGACAAAAAGAAGCCAGTTGTTGTTTCAGGTCCATAT
ATAGTTTATATGCTACAAGATCTTGATATTCTGAAGACTGGACAACAATTAGGAAGGCAATGGCTACAT
TGGGGCCACACAGAGTGAACGGAACCACTGTGAACTGGAAAAACATCTGCACAGTGCTAGATCTGA
AGAGGGAAGACTATATTATGATGGTGAATGGTATATACGTGGACAAACAATATGTATTGATAAAAAAGAT
GAATGTCCTACAAGTGTGTAATTACAACAATTAACCATGATGAAGTTTGGTTTAAGAGGCCTGATGGAA
GCAATCTAAGCTTTACATTTACAGCTACAGAAAGGAAAATATTCAATTAACATTCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202645 protein sequence Red=Cloning site Green=Tags(s) MPVHSRGDKKETNHHDEMEVDYAENEGSSSEDEDTESSSVSEGDGSSEMDDDCERRRMECLDEMSNLEK QFTDLKDQLYKERLSQVDAKLQEVIAGKAPEYLEPLATLQENMQIRTKVAGIYRELCLESVKNKYECEIQ ASRQHCESEKLLLYDTVQSELEEKIRRLEEDRHSIDITSELWNDELQSRKKRKDPFSPDKKKPVVSGPY IVYMLQDLILEDWTTIRKAMATLGPHRVKTEPPVKLEKHLHSARSEEGRLYYDGEWYIRGQTICIDKKD ECPTSAVITTTINHDEVWFKRPDGSKSKLYISQLQKGKYSIKHS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_032352
ORF Size:	969 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032352.3</u>
RefSeq Size:	2604 bp
RefSeq ORF:	972 bp
Locus ID:	84312
UniProt ID:	<u>Q5PSV4</u>
Cytogenetics:	14q13.2
MW:	37.6 kDa