

## Product datasheet for **RC202652L3V**

### Triosephosphate isomerase (TPI1) (NM\_000365) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Triosephosphate isomerase (TPI1) (NM_000365) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Triosephosphate isomerase
Synonyms:	HEL-S-49; TIM; TPI; TPID
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000365
ORF Size:	747 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202652).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000365.4</a>
RefSeq Size:	1366 bp
RefSeq ORF:	750 bp
Locus ID:	7167
UniProt ID:	<a href="#">P60174</a>
Cytogenetics:	12p13.31
Domains:	TIM



[View online »](#)

<b>Protein Pathways:</b>	Fructose and mannose metabolism, Glycolysis / Gluconeogenesis, Inositol phosphate metabolism, Metabolic pathways
<b>MW:</b>	26.7 kDa
<b>Gene Summary:</b>	This gene encodes an enzyme, consisting of two identical proteins, which catalyzes the isomerization of glyceraldehydes 3-phosphate (G3P) and dihydroxy-acetone phosphate (DHAP) in glycolysis and gluconeogenesis. Mutations in this gene are associated with triosephosphate isomerase deficiency. Pseudogenes have been identified on chromosomes 1, 4, 6 and 7. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2009]