

## Product datasheet for **RC202662L4V**

### THYN1 (NM\_001037304) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	THYN1 (NM_001037304) Human Tagged ORF Clone Lentiviral Particle
Symbol:	THYN1
Synonyms:	HSPC144; MDS012; MY105; THY28; THY28KD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001037304
ORF Size:	498 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202662).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001037304.1</a>
RefSeq Size:	790 bp
RefSeq ORF:	501 bp
Locus ID:	29087
UniProt ID:	<a href="#">Q9P016</a>
Cytogenetics:	11q25
MW:	18.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that is highly conserved among vertebrates and plant species and may be involved in the induction of apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]