

Product datasheet for **RC202672L4V**

METTL2A (NM_001005372) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	METTL2A (NM_001005372) Human Tagged ORF Clone Lentiviral Particle
Symbol:	METTL2A
Synonyms:	FLJ12760; methyltransferase like 2A; METTL2; METTL2, FLJ12760
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001005372
ORF Size:	726 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202672).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001005372.1 , NP_001005372.1
RefSeq Size:	2705 bp
RefSeq ORF:	728 bp
Locus ID:	339175
Cytogenetics:	17q23.2
Protein Families:	Druggable Genome
MW:	27.5 kDa



[View online »](#)

Gene Summary:

Probable S-adenosyl-L-methionine-dependent methyltransferase that mediates 3-methylcytidine modification of some tRNAs.[UniProtKB/Swiss-Prot Function]