

Product datasheet for **RC202691L4V**

EIF1B (NM_005875) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EIF1B (NM_005875) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EIF1B
Synonyms:	GC20
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005875
ORF Size:	339 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202691).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005875.2
RefSeq Size:	1014 bp
RefSeq ORF:	342 bp
Locus ID:	10289
UniProt ID:	O60739
Cytogenetics:	3p22.1
Domains:	SUI1
MW:	12.8 kDa



[View online »](#)

Gene Summary: Probably involved in translation.[UniProtKB/Swiss-Prot Function]