

Product datasheet for **RC202692L4V**

G3BP (G3BP1) (NM_005754) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	G3BP (G3BP1) (NM_005754) Human Tagged ORF Clone Lentiviral Particle
Symbol:	G3BP
Synonyms:	G3BP; HDH-VIII
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005754
ORF Size:	1398 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202692).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005754.2
RefSeq Size:	2855 bp
RefSeq ORF:	1401 bp
Locus ID:	10146
UniProt ID:	Q13283
Cytogenetics:	5q33.1
Domains:	RRM, NTF2
Protein Families:	Druggable Genome




[View online »](#)

MW: 52.2 kDa

Gene Summary: This gene encodes one of the DNA-unwinding enzymes which prefers partially unwound 3'-tailed substrates and can also unwind partial RNA/DNA and RNA/RNA duplexes in an ATP-dependent fashion. This enzyme is a member of the heterogeneous nuclear RNA-binding proteins and is also an element of the Ras signal transduction pathway. It binds specifically to the Ras-GTPase-activating protein by associating with its SH3 domain. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Jul 2008]

Product images:



[RC202692L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC202692L4V particle to overexpress human G3BP1-mGFP fusion protein.