

Product datasheet for **RC202696A1V**

Human CRSP9 (MED7) (NM_004270) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CRSP9 (MED7) (NM_004270) AAV Particle
Tag: Myc-DDK
Symbol: CRSP9
Synonyms: ARC34; CRSP9; CRSP33
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202696 representing NM_004270
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGTGAACACAGCAAGTGAGTGCACCTCCACCACCTCCAATGCAATATATCAAGGAATATACGGATG
AAAATATTCAAGAAGGCTTAGCTCCCAAGCCTCCCCCTCCAATAAAAGACAGTTACATGATGTTTGGCAA
TCAGTTCCAATGTGATGATCTTATCATCCGCCCTTTGGAAAGTCAGGGCATCGAACGGCTTCATCCTATG
CAGTTTGATCACAAGAAAGAACTGAGAAAATTAATATGTCTATCCTTATTAATTTCTTGACCTTTTAG
ATATTTAATAAGGAGCCCTGGGAGTATAAACGAGAAGAGAACTAGAAGATCTTAAGCTGCTTTTGT
ACACGTGCATCATCTTATAAATGAATACCGACCCACCAAGCAAGAGAGACCTTGAGAGTCATGATGGAG
GTCCAGAAACGTCAACGGCTTGAACAGCTGAGAGATTTCAAAAGCACCTGGAACGAGTAATTGAAATGA
TTCAGAAATGCTTGGCTTCTTTGCCTGATGATTTGCCTCATTGAGAGCAGGAATGAGAGTAAAACTGA
ACCAATGGATGCTGATGATAGCAACAATTGACTGGACAGAATGAACATCAAAGAGAAAATTCAGGTCAT
AGGAGAGATCAGATTATAGAGAAAGATGCTGCCTTGTGTCTTAATTGATGAGATGAATGAAAGACCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202696 representing NM_004270
Red=Cloning site Green=Tags(s)

MGEPPQVSALPPPPMQYIKEYTDENIQEGLAPKPPPIKDSYMMFGNQFQCDDLIIIRPLESQGIERLHPM
QFDHKKELRKLNMSILINFLDLLDILIRSPGSIKREEKLEDLKLFFVHVHHLINERYPHQARETLRVMM
VQKRQRLETAERFQKHLERVIEMIQNCLASLPDDLPHSEAGMRVKTEPMDADDSNNCTGQNEHQRENSGH
RRDQIIIEKDAALCVLIDEMNERP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_004270
ORF Size:	699 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004270.3
RefSeq Size:	1066 bp
RefSeq ORF:	702 bp
Locus ID:	9443
UniProt ID:	O43513
Cytogenetics:	5q33.3
MW:	27.1 kDa