

## Product datasheet for **RC202699L1V**

### GGPS1 (NM\_004837) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | Lentiviral Particles                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Product Name:             | GGPS1 (NM_004837) Human Tagged ORF Clone Lentiviral Particle                                                                                                                                                                                                                                                                                                                                                                                       |
| Symbol:                   | GGPS1                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Synonyms:                 | GGPPS; GGPPS1                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Mammalian Cell Selection: | None                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Vector:                   | pLenti-C-Myc-DDK (PS100064)                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Tag:                      | Myc-DDK                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ACCN:                     | NM_004837                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| ORF Size:                 | 900 bp                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC202699).                                                                                                                                                                                                                                                                                                                                                                                     |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| RefSeq:                   | <a href="#">NM_004837.3</a>                                                                                                                                                                                                                                                                                                                                                                                                                        |
| RefSeq Size:              | 2921 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| RefSeq ORF:               | 903 bp                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Locus ID:                 | 9453                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Cytogenetics:             | 1q42.3                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Domains:                  | polyprenyl_synt                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Protein Pathways:         | Metabolic pathways, Terpenoid backbone biosynthesis                                                                                                                                                                                                                                                                                                                                                                                                |
| MW:                       | 34.9 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                           |



[View online »](#)

**Gene Summary:**

This gene is a member of the prenyltransferase family and encodes a protein with geranylgeranyl diphosphate (GGPP) synthase activity. The enzyme catalyzes the synthesis of GGPP from farnesyl diphosphate and isopentenyl diphosphate. GGPP is an important molecule responsible for the C20-prenylation of proteins and for the regulation of a nuclear hormone receptor. Alternate transcriptional splice variants, both protein-coding and non-protein-coding, have been found for this gene. [provided by RefSeq, Sep 2010]