

Product datasheet for **RC202715L4V**

NDUFA2 (NM_002488) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NDUFA2 (NM_002488) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NDUFA2
Synonyms:	B8; CD14; CIB8; MC1DN13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002488
ORF Size:	297 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202715).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002488.2
RefSeq Size:	726 bp
RefSeq ORF:	300 bp
Locus ID:	4695
UniProt ID:	O43678
Cytogenetics:	5q31.3
Domains:	L51_S25_CI-B8



[View online »](#)

Protein Pathways:	Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease
MW:	10.9 kDa
Gene Summary:	The encoded protein is a subunit of the hydrophobic protein fraction of the NADH:ubiquinone oxidoreductase (complex 1), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane, and may be involved in regulating complex I activity or its assembly via assistance in redox processes. Mutations in this gene are associated with Leigh syndrome, an early-onset progressive neurodegenerative disorder. Alternative splicing results in multiple transcript variants.[provided by RefSeq, May 2010]