

Product datasheet for **RC202729A1V**

Human TNFAIP8 (NM_014350) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TNFAIP8 (NM_014350) AAV Particle
Tag: Myc-DDK
Symbol: TNFAIP8
Synonyms: GG2-1; MDC-3.13; NDED; SCC-S2; SCCS2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202729 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCACTCCGAAGCAGAAGAATCCAAGGAAGTGGCCACAGATGTCTTTAATTCCAAAACCTGGCCGTTCC
AGGCACAAAAGAAGATCTTGGGTAAATGGTGTCCAATCCATCGCCACCACCTTAATAGACGACACAAG
TAGTGAGGTGCTGGATGAGCTCTACAGAGTGACCAGGGAGTACACCCAAAACAAGAAGGAGGCAGAGAAG
ATCATCAAGAACCTCATCAAGACAGTCATCAAGCTGGCCATTCTTTATAGGAATAATCAGTTTAATCAAG
ATGAGCTAGCATTGATGGAGAAATTTAAGAAGAAAGTTCATCAGCTTGCTATGACCGTGGTCAGTTTCCA
TCAGGTGATTATACCTTTGACCGGAATGTGTTATCCAGGCTGTAAATGAATGCAGAGAGATGCTGCAC
CAAATCATTACAGCGCCACCTCACTGCCAAGTCACATGGACGGGTAATAATGTGTTTGATCATTTTTTCAG
ATTGTGAATTTTTGGCTGCCTTGATAATCCTTTTGGGAATTTTAAACCCCACTTACAAAACTATGTGA
TGGTATCAACAAATGTTGGATGAAGAGAACATA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202729 protein sequence
Red=Cloning site Green=Tags(s)

MHSEAEESKEVATDVFN SKNLAVQAQKKILGKMVSKSIATTLIDDTSSSEVLDEL YRV TREYTQNKKEAEK
IIKNLIKTVIKLAILYRNNQFNQDELALMEKFKKKVHQLAMTVVSFHHQVDYTFDRNVL SLLNECREMLH
QIIQRHLTAKSHGRVNNVDFHFSDFLAALYNPFGNFKPHLQKLCDGINKMLDEENI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_014350
ORF Size:	594 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_014350.2
RefSeq Size:	2103 bp
RefSeq ORF:	597 bp
Locus ID:	25816
UniProt ID:	O95379
Cytogenetics:	5q23.1
MW:	23 kDa