

## Product datasheet for **RC202739L1V**

### **FXYD2 (NM\_021603) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	FXYD2 (NM_021603) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FXYD2
Synonyms:	ATP1G1; HOMG2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_021603
ORF Size:	192 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202739).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021603.2</a>
RefSeq Size:	591 bp
RefSeq ORF:	195 bp
Locus ID:	486
UniProt ID:	<a href="#">P54710</a>
Cytogenetics:	11q23.3
Protein Families:	Druggable Genome, Ion Channels: Other, Transmembrane
MW:	7.4 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the FXVD family of transmembrane proteins. This particular protein encodes the sodium/potassium-transporting ATPase subunit gamma. Mutations in this gene have been associated with Renal Hypomagnesemia-2. Alternatively spliced transcript variants have been described. Read-through transcripts have been observed between this locus and the upstream FXVD domain-containing ion transport regulator 6 (FXVD6, GeneID 53826) locus.[provided by RefSeq, Feb 2011]