

## Product datasheet for **RC202756L4V**

### CD69 (NM\_001781) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CD69 (NM_001781) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CD69
Synonyms:	AIM; BL-AC/P26; CLEC2C; EA1; GP32/28; MLR-3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001781
ORF Size:	597 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202756).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001781.1</a>
RefSeq Size:	1676 bp
RefSeq ORF:	600 bp
Locus ID:	969
UniProt ID:	<a href="#">Q07108</a>
Cytogenetics:	12p13.31
Domains:	CLECT
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane


[View online »](#)

**MW:** 22.6 kDa

**Gene Summary:** This gene encodes a member of the calcium dependent lectin superfamily of type II transmembrane receptors. Expression of the encoded protein is induced upon activation of T lymphocytes, and may play a role in proliferation. Furthermore, the protein may act to transmit signals in natural killer cells and platelets. [provided by RefSeq, Aug 2011]