

Product datasheet for **RC202778L3V**

FGL1 (NM_201553) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | FGL1 (NM_201553) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | FGL1 |
| Synonyms: | HFREP1; HP-041; HPS; LFIRE-1; LFIRE1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_201553 |
| ORF Size: | 936 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC202778). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_201553.1 , NP_963847.1 |
| RefSeq Size: | 1500 bp |
| RefSeq ORF: | 939 bp |
| Locus ID: | 2267 |
| UniProt ID: | Q08830 |
| Cytogenetics: | 8p22 |
| Protein Families: | Druggable Genome, Secreted Protein |
| MW: | 36.4 kDa |



[View online »](#)

Gene Summary:

Fibrinogen-like 1 is a member of the fibrinogen family. This protein is homologous to the carboxy terminus of the fibrinogen beta- and gamma- subunits which contains the four conserved cysteines of fibrinogens and fibrinogen related proteins. However, this protein lacks the platelet-binding site, cross-linking region and a thrombin-sensitive site which are necessary for fibrin clot formation. This protein may play a role in the development of hepatocellular carcinomas. Four alternatively spliced transcript variants encoding the same protein exist for this gene. [provided by RefSeq, Jul 2008]