

## **Product datasheet for RC202786L1**

PSMA4 (NM\_002789) Human Tagged Lenti ORF Clone

### 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsusport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

# Product data:

**Product Type:** Expression Plasmids

**Product Name:** PSMA4 (NM\_002789) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: PSMA4

Synonyms: HC9; HsT17706; PSC9

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC202786).

Sequence:

**Restriction Sites:** Sgfl-Mlul

Cloning Scheme:





<sup>\*</sup> The last codon before the Stop codon of the ORF.

**ACCN:** NM\_002789

ORF Size: 783 bp





#### PSMA4 (NM\_002789) Human Tagged Lenti ORF Clone - RC202786L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 002789.3</u>

 RefSeq Size:
 1235 bp

 RefSeq ORF:
 786 bp

 Locus ID:
 5685

 UniProt ID:
 P25789

Cytogenetics: 15q25.1

**Domains:** proteasome

**Protein Families:** Druggable Genome, Protease


Protein Pathways: Proteasome MW: 29.5 kDa

**Gene Summary:** This gene encodes a core alpha subunit of the 20S proteosome, which is a highly ordered


ring-shaped structure composed of four rings of 28 non-identical subunits. Proteasomes cleave peptides in an ATP- and ubiquitin-dependent manner. [provided by RefSeq, Aug 2016]



## **Product images:**



Circular map for RC202786L1



Double digestion of RC202786L1 using Sgfl and Mlul  $\,$