

Product datasheet for **RC202787A1V**

Human ARHI (DIRAS3) (NM_004675) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human ARHI (DIRAS3) (NM_004675) AAV Particle
Tag: Myc-DDK
Symbol: ARHI
Synonyms: ARHI; NOEY2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202787 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGTAAAGCCAGCTTTGGCTCCAAGGAACAGAAGCTGCTGAAGCGGTTGCGGCTTCTGCCCGCCCTGC
TTATCCTCCGCGCCTTCAAGCCCCACAGGAAGATCAGAGATTACCGCGTCGTGGTAGTCGGCACCCTGG
TGTGGGAAAAAGTACGCTGCTGCACAAGTGGGCGAGCGGCAACTTCCGTCATGAGTACCTGCCGACCAT
GAAAATACCTACTGCCAGTTGCTGGGCTGCAGCCACGGTGTGCTTCCCTGCACATCACCAGACGAAGA
GTGGCGACGGCAACCGCGCTCTGCAGCGCCACGTTATAGCCCGGGGCCACGCCTTCGTCCTGGTCTACTC
AGTCACCAAGAAGGAAACCTTGAAGAGCTGAAGGCCTTCTATGAGCTGATCTGCAAGATCAAAGGTAAC
AACCTGCATAAGTTCCCATCGTGTGGTGGGCAATAAAGTGATGACACCCACCGGGAGGTGGCCTGA
ATGATGGTGCCACCTGTGCGATGGAGTGAATTGCGCCTTCATGGAGATTTAGCCAAGACCGATGTGAA
TGTGCAGGAGCTGTTCCACATGCTGCTGAATTACAAGAAAAAGCCACACCGGCCTCCAGGAGCCCGAG
AAGAAATCCCAGATGCCAACACCACTGAGAAGCTGCTTGACAAGTGCATAATCATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202787 protein sequence
Red=Cloning site Green=Tags(s)

MGNASFGSKEQKLLKRLRLPALILRAFKPHRKIRDYRVVVVGTAGVGKSTLLHKWASGNFRHEYLP
ENTYQQLLGCSHGVL SLHITDSKSGDGNRALQRHVIARGHAFVLVYSVTKKETLEELKAFYELICKIKGN
NLHKFPIVLVGNKSDDTHREVALNDGATCAMEWNAFMEISAKTDVNVQELFHMLLNYKKKPTTGLQEPE
KKSQMPNTTEKLLDKCIIM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_004675
ORF Size:	687 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004675.2
RefSeq Size:	1642 bp
RefSeq ORF:	690 bp
Locus ID:	9077
UniProt ID:	O95661
Cytogenetics:	1p31.3
MW:	25.9 kDa