

Product datasheet for **RC202793L3V**

LSM3 (NM_014463) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | LSM3 (NM_014463) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | LSM3 |
| Synonyms: | SMX4; USS2; YLR438C |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_014463 |
| ORF Size: | 306 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC202793). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_014463.1 |
| RefSeq Size: | 695 bp |
| RefSeq ORF: | 309 bp |
| Locus ID: | 27258 |
| UniProt ID: | P62310 |
| Cytogenetics: | 3p25.1 |
| Domains: | Sm |
| Protein Families: | Stem cell - Pluripotency |



[View online »](#)

Protein Pathways: RNA degradation, Spliceosome

MW: 11.8 kDa

Gene Summary: Sm-like proteins were identified in a variety of organisms based on sequence homology with the Sm protein family (see SNRPD2; MIM 601061). Sm-like proteins contain the Sm sequence motif, which consists of 2 regions separated by a linker of variable length that folds as a loop. The Sm-like proteins are thought to form a stable heteromer present in tri-snRNP particles, which are important for pre-mRNA splicing.[supplied by OMIM, Apr 2004]