

Product datasheet for **RC202809A1V**

Human PSMC6 (NM_002806) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PSMC6 (NM_002806) AAV Particle
Tag:	Myc-DDK
Symbol:	PSMC6
Synonyms:	p42; RPT5; SUG2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202809 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGACCCTAGAGATAAGGCGCTTCAGGACTACCGCAAGAAGTTGCTTGAACACAAGGAGATCGACG
GCCGTCTTAAGGAGTTAAGGGAACAATTAAGAAGCTTACCAAGCAGTATGAAAAGTCTGAAAATGATCT
GAAGGCCCTACAGAGTGTGGGCAGATCGTGGGTGAAGTGCTTAACAGTTAACTGAAGAAAAATTCATT
GTTAAAGCTACCAATGGACCAAGATATGTTGTGGTTGTCGTCGACAGCTTGACAAAAGTAAGCTGAAGC
CAGGAACAAGAGTTGCTTTGGATATGACTACCTAACTATCATGAGATATTTGCCGAGAGAGGTGGATCC
ACTGGTTTATAACATGTCTCATGAGGACCCTGGGAATGTTTCTTATTCTGAGATTGGAGGGCTATCAGAA
CAGATCCGGGAATTAAGAGAGGTGATAGAATTACCTCTTACAAACCCAGAGTTATTTAGCGGTGATAGGAA
TAATACCTCCAAAAGGCTGTTTGTATATGGACCACCAGGTACGGGAAAAACACTCTTGGCACGAGCCGT
TGCTAGCCAGCTGGACTGCAATTTCTAAAGGTTGTATCTAGTTCTATTGTAGACAAGTACATTGGTGAA
AGTGCTCGTTTGATCAGAGAAATGTTAATTATGCTAGAGATCATCAACCATGCATCATTTTATGGATG
AAATAGATGCTATTGGTGGTCGTCGTTTCTGAGGGTACTTCAGCTGACAGAGAGATTCAGAGAACGTT
AATGGAGTTACTGAATCAAATGGATGGATTTGATACTCTGCATAGAGTTAAATGATCATGGCTACAAAC
AGACCAGATACACTGGATCCTGCTTTGCTGCGTCCAGGAAGATTAGATAGAAAAATACATATTGATTTGC
CAAATGAACAAGCAAGATTAGACATACTGAAAATCCATGCAGGTCCCATTACAAAGCATGGTGAAATAGA
TTATGAAGCAATTGTGAAGCTTTCGGATGGCTTAAATGGAGCAGATCTGAGAAATGTTTGTACTGAAGCA
GGTATGTTTCGCAATTCGTGCTGATCATGATTTTGTAGTACAGGAAGACTTCATGAAAGCAGTCAGAAAAG
TGGCTGATTCTAAGAAGCTGGAGTCTAAATTGGACTACAACTGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC202809 protein sequence Red=Cloning site Green=Tags(s) MADPRDKALQDYRKLLLEHKEIDGRLKELREQLKELTKQYEKSENDLKALQSVGQIVGEVLKQLTEEKFI VKATNGPRYVVGCRRLDKSKLKPGRVALDMTTLTIMRYLPREVDPLVYNMSHEDPGNVSYSYSEIGGLSE QIRELREVIPLTNPELFQRVGIIPKGCCLLYGPPGTGKTLLARAVASQLDCNFLKVVSSSIVDKYIGE SARLIREFNYARDHQPCIIIFMDEIDAIGGRRFSEGTSADREIQRTLMELLNQMDGFDLHRVKMIMATN RPDTLDPALLRPGRLDRKIHIDLPNEQARLDILKIHAGPITKHGEIDYEAIVKLSDFNGADLRNVCTEA GMFAIRADHDFVVQEDFMKAVRKVADSKKLESKLDYKPV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002806
ORF Size:	1167 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_002806.2
RefSeq Size:	1599 bp
RefSeq ORF:	1170 bp
Locus ID:	5706
UniProt ID:	P62333
Cytogenetics:	14q22.1
MW:	44.2 kDa