

Product datasheet for **RC202817A1V**

Human DPM1 (NM_003859) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DPM1 (NM_003859) AAV Particle
Tag:	Myc-DDK
Symbol:	DPM1
Synonyms:	CDGIE; MPDS
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202817 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCTCCTTGAAGTCAGTCGTAGTCCTCGCAGGTCTCGGCGGGAGCTGGAAGTGCAGTCCACGAC
AGAACAAATATTCGGTGCTTTTACCTACCTACAACGAGCGCGAGAACCTGCCGCTCATCGTGTGGCTGCT
GGTGAAAAGCTTCTCCGAGAGTGGAATCAACTATGAAATTATAATCATAGATGATGGAAGCCAGATGGA
ACAAGGGATGTTGCTGAACAGTTGGAGAAGATCTATGGGTGAGACAGAATTCTTCTAAGACCACGAGAGA
AAAAGTTGGGACTAGGAATGCATATATTCATGGAATGAAACATGCCACAGGAACTACATCATTATTAT
GGATGCTGATCTCTCACACCATCCAAAATTTATTCCTGAATTTATTAGGAAGCAAAAGGAGGGTAATTTT
GATATTGTCTCTGGAATCGCTACAAAGGAAATGGAGGTGTATATGGCTGGGATTGAAAAGAAAATAA
TCAGCCGTGGGGCCAATTTTTTAACTCAGATCTTGCTGAGACCAGGAGCATCTGATTTAACAGGAAGTTT
CAGATTATACGAAAAGAAGTTCTAGAGAAATTAATAGAAAAATGTGTTTCTAAAGGCTACGTCTTCCAG
ATGGAGATGATTGTTCCGGCAAGACAGTTGAATTATACTATTGGCGAGGTTCCAATATCATTTGTGGATC
GTGTTTATGGTGAATCCAAGTTGGGAGGAAATGAAATAGTATCTTTCTTGAAAGGATTATTGACTCTTTT
TGCTACTACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202817 protein sequence Red=Cloning site Green=Tags(s) MASLEVSRSRRSRRELEVRSPRQNKYSVLLPTYNERENLPLIVLLVKSFSESGINYEIIIIDDGSPDG TRDVAEQLEKIYGSDRILLRPREKKLGLGTAYIHGMKHATGNYIIIMDADLSHHPKFIPEFIRKQKEGNF DIVSGTRYKGNNGVYGDWLKRKIIISRGANFLTQILLRPGASDLTGSFRLYRKEVLEKLEKCVSKGYVFQ MEMIVRARQLNYTIGEVPISFVDRVYGESKLGNEIVSFLKGLLTFLFATT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_003859
ORF Size:	780 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003859.1</u>
RefSeq Size:	1047 bp
RefSeq ORF:	783 bp
Locus ID:	8813
UniProt ID:	<u>O60762</u>
Cytogenetics:	20q13.13
MW:	29.6 kDa