

Product datasheet for **RC202820L4V**

C11orf2 (VPS51) (NM_013265) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C11orf2 (VPS51) (NM_013265) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C11orf2
Synonyms:	ANG2; ANG3; C11orf2; C11orf3; FFR; PCH13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013265
ORF Size:	2346 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202820).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013265.2
RefSeq Size:	2546 bp
RefSeq ORF:	2349 bp
Locus ID:	738
UniProt ID:	Q9UID3
Cytogenetics:	11q13.1
MW:	85.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the vacuolar protein sorting-associated protein 51 family. The encoded protein is a component of the Golgi-associated retrograde protein complex which acts as a tethering factor for carriers in retrograde transport from the early and late endosomes to the trans-Golgi network. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012]