

Product datasheet for **RC202825A1V**

Human CIAO2B (NM_016062) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CIAO2B (NM_016062) AAV Particle
Tag: Myc-DDK
Symbol: CIAO2B
Synonyms: CGI-128; CIA2B; FAM96B; MIP18
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202825 ORF sequence
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTAGGCGGCGGGGGTTCGGCGGCGGCTCCTGGAGAATGCCAACCCCTCATCTACCAGCGCTCTG
GGGAGCGGCTGTGACGGCAGGCGAGGAGGACGAGCAGGTTCCCGACAGCATCGACGCACGCGAGATCTT
CGATCTGATTGCTCCATCAATGACCCGGAGCATCCACTGACGCTAGAGGAGTTGAACGTAGTAGAGCAG
GTGCGGGTTTCAGGTTAGCGACCCCGAGAGTACAGTGGCTGTGGCTTTACACCAACCATTCCGCACTGCA
GCATGGCCACCCTTATTGGTCTGTCCATCAAGTCAAGCTTCTGCGCTCCCTTCCTCAGCGTTTCAAGAT
GGACGTGCACATTACTCCGGGGACCCATGCCTCAGAGCATGCAGTGAACAAGCAACTGCAGATAAGGAG
CGGGTGGCAGCTGCCCTGGAGAACCCACCTCTTGGAGGTTGTGAATCAGTGCCTGTCAGCCCGCTCC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202825 protein sequence
Red=Cloning site **Green**=Tags(s)

MYGGGVGGGLENANPLIYQSGERPVTAGEEDEQVPDSIDAREIFDLIRSINDPEHPLTLEELNVVEQ
VRVQVSDPESTVAVFTPTIPHCSMATLIGLSIKVKLLRSLPQRFKMDVHITPGTHASEHAVNKQLADKE
RVAAALENTHLLEVNNQCLSARS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_016062
ORF Size:	489 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_016062.1</u>
RefSeq Size:	716 bp
RefSeq ORF:	492 bp
Locus ID:	51647
UniProt ID:	<u>Q9Y3D0</u>
Cytogenetics:	16q22.1
MW:	17.7 kDa