

Product datasheet for **RC202826A1V**

Human UBC6e (UBE2J1) (NM_016021) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human UBC6e (UBE2J1) (NM_016021) AAV Particle
Tag:	Myc-DDK
Symbol:	UBC6e
Synonyms:	CGI-76; HSPC153; HSPC205; HSU93243; NCUBE-1; NCUBE1; UBC6; UBC6E; Ubc6p
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202826 representing NM_016021 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGACCCGCTACAACCTGAAGAGTCCGGCTGTAAACGTTTAAATGAAAGAAGCGGCAGAATTGAAAG
ATCCAACAGATCATTACCATGCGCAGCCTTTAGAGGATAACCTTTTGAATGGCACTTCACGGTTAGAGG
GCCCCAGACTCCGATTTTGATGGAGGAGTTTATCACGGGCGGATAGTACTGCCACCAGAGTATCCCATG
AAACCACCAAGCATTATTCTCTAACGGCTAATGGTCGATTTGAAGTGGGCAAGAAAATCTGTTTGAGCA
TCTCAGGCCATCATCTGAAACTTGGCAGCCTTCGTGGAGTATAAGGACAGCATTATTAGCCATCATTGG
GTTTATGCCAACAAGGAGAGGGAGCCATAGGTTCTCTAGATTACACTCCTGAGGAAAGAAGAGCACTT
GCCAAAAAATCACAAGATTTCTGTTGTGAAGGATGTGGCTCTGCCATGAAGGATGTCTGTTGCCTTTAA
AATCTGGAAGCGATTCAAGCCAAGCTGACCAAGAAGCCAAGAACTGGCTAGGCAAAATAAGCTTTAAGGC
AGAAGTCAATTCATCTGGAAAGACTATCTCTGAGTCAGACTTAAACCACTCTTTTCACTAACTGATTTA
CAAGATGATATACCTACAACATTCCAGGGTGCTACGGCCAGTACATCGTACGGAGTCCAGAATTCCTCAG
CAGCATCCTTTCATCAACCTACCAACCTGTAGCTAAGAATACCTCCATGAGCCCTCGACAGCGCCGGGC
CCAGCAGCAGAGTCAGAGAAGGTTGTCTACTTCACCAGATGTAATCCAGGGCCACCAGCCAAGAGACAAC
CACACTGATCATGGTGGGTGAGCTGTACTGATTGTCTCCTGACTTTGGCATTGGCAGCTCTTATATTCC
GACGAATATATCTGGCAAACGAATACATATTTGACTTTGAGTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202826 representing NM_016021 Red=Cloning site Green=Tags(s) METRYNLKSPAVKRLMKEAAELKDPTDHYHAQPLEDNLFWEHFTVRGPPDSDFDGGVYHGRIVLPPEYPM KPPSIILLTANGRFEVGKKICLSISGHPETWQPSWSIRTALLAIIGFMPTKGEGAIGSLDYTPERRAL AKKSQDFCCEGCGSAMKDVLLPLKSGSDSSQADQEAKELARQISFKAENVSSGKTISESDLNHSFSLTDL QDDIPTTFQGATASTSYGVQNSSAASFHQPTQPVAKNTSMSPRQRRQQSQRRRLSTSPDVIQGHQPRDN HTDHGGS AVLIVILTLALAALIFRRIYLANEYIFDFEL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_016021
ORF Size:	954 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_016021.2 , NP_057105.2
RefSeq Size:	4360 bp
RefSeq ORF:	957 bp
Locus ID:	51465
UniProt ID:	Q9Y385
Cytogenetics:	6q15
MW:	35 kDa