

Product datasheet for **RC202840L4V**

NUDT11 (NM_018159) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NUDT11 (NM_018159) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NUDT11
Synonyms:	APS1; ASP1; DIPP3b; DIPP3beta; hDIPP3beta
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018159
ORF Size:	492 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202840).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018159.2
RefSeq Size:	2413 bp
RefSeq ORF:	495 bp
Locus ID:	55190
UniProt ID:	Q96G61
Cytogenetics:	Xp11.22
Domains:	NUDIX
MW:	18.6 kDa



[View online »](#)

Gene Summary:

NUDT11 belongs to a subgroup of phosphohydrolases that preferentially attack diphosphoinositol polyphosphates (Hidaka et al., 2002 [PubMed 12105228]).[supplied by OMIM, Mar 2008]