

## Product datasheet for **RC202850L2V**

### **NENF (NM\_013349) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	NENF (NM_013349) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NENF
Synonyms:	CIR2; SCIRP10; SPUF
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_013349
ORF Size:	516 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202850).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013349.3</a>
RefSeq Size:	944 bp
RefSeq ORF:	519 bp
Locus ID:	29937
UniProt ID:	<a href="#">Q9UMX5</a>
Cytogenetics:	1q32.3
Domains:	heme_1
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

**MW:** 18.86 kDa

**Gene Summary:** This gene encodes a neurotrophic factor that may play a role in neuron differentiation and development. A pseudogene of this gene is found on chromosome 12. Alternate splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jan 2009]