

Product datasheet for **RC202851L4V**

Hey L (HEYL) (NM_014571) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hey L (HEYL) (NM_014571) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Hey L
Synonyms:	bHLHb33; HESR3; HEY3; HRT3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014571
ORF Size:	984 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202851).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014571.2
RefSeq Size:	4131 bp
RefSeq ORF:	987 bp
Locus ID:	26508
UniProt ID:	Q9NQ87
Cytogenetics:	1p34.2
Protein Families:	Druggable Genome, Transcription Factors
MW:	35.1 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the hairy and enhancer of split-related (HESR) family of basic helix-loop-helix (bHLH)-type transcription factors. The sequence of the encoded protein contains a conserved bHLH and orange domain, but its YRPW motif has diverged from other HESR family members. It is thought to be an effector of Notch signaling and a regulator of cell fate decisions. Alternatively spliced transcript variants have been found, but their biological validity has not been determined. [provided by RefSeq, Jul 2008]