

Product datasheet for **RC202852L3V**

GARS (GARS1) (NM_002047) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GARS (GARS1) (NM_002047) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GARS1
Synonyms:	CMT2D; DSMAV; GARS; GlyRS; HMN5; HMN5A; SMAD1; SMAJI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002047
ORF Size:	2217 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202852).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002047.2
RefSeq Size:	2759 bp
RefSeq ORF:	2220 bp
Locus ID:	2617
UniProt ID:	P41250
Cytogenetics:	7p14.3
Domains:	WHEP-TRS, tRNA-synt_2b, HGTP_anticodon
Protein Pathways:	Aminoacyl-tRNA biosynthesis



[View online »](#)

MW: 83.1 kDa

Gene Summary: This gene encodes glycyl-tRNA synthetase, one of the aminoacyl-tRNA synthetases that charge tRNAs with their cognate amino acids. The encoded enzyme is an (alpha)₂ dimer which belongs to the class II family of tRNA synthetases. It has been shown to be a target of autoantibodies in the human autoimmune diseases, polymyositis or dermatomyositis. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2015]