

Product datasheet for **RC202856A1V**

Human ZNF346 (NM_012279) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ZNF346 (NM_012279) AAV Particle
Tag:	Myc-DDK
Symbol:	ZNF346
Synonyms:	JAZ; Zfp346
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202856 representing NM_012279 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGTATCCCGCGCCGGCCACGGTGCAGGCCGCGGACGGCGGAGCGGCCGGGCCTTACAGCAGCTCGG
AGTTGCTGGAGGGCCAGGAGCCGGACGGGTGCGCTTTGACCGCGAGAGGGCGCGCCGCTGTGGGAAGC
CGTGTCGGTGCCCGAGCCGGTGGGTAGAGAGGAAGTGGAGCACATGATCCAGAAGAACCAATGTCTCTTC
ACCAACACCCAGTGTAAGGTTTGCTGCGCCTTGCTTATTTCTGAGTCCCAGAAGCTGGCACATTACCAGA
GCAAAAAACATGCCAACAAAGTGAAGAGATACCTAGCAATCCATGGAATGGAGACATTAAAGGGGAAAC
GAAGAAGCTAGACTCAGATCAGAAGAGCAGCAGAAGCAAGACAAGAACCAGTGCTGCCCATCTGTAAAC
ATGACCTTTTCTCCCTGTCGTGGCCAGTCGCACTACCTGGGGAAGACCCACGCAAGAACTTAAAGC
TGAAGCAGCAGTCCACTAAGGTGGAAGCCTTGACCAGAATAGAGAGATGATAGACCCAGACAAGTTCTG
CAGCCTCTGCCATGCACTTTCAACGACCCTGTCATGGCTCAACAACATTATGTGGGCAAGAAACACAGA
AAACAGGAGACCAAGCTCAAACATATGGCACGCTATGGCGGCTGGCGGACCCTGCTGCTACTGACTTTC
CAGCTGGAAGGGCTACCCCTGCAAAACATGTAAGATAGTGCTGAATCCATAGAACAGTACCAAGCTCA
TGTCAGCGGCTTCAAACACAAGAACCAGTCACCAAAAACAGTGGCATCATCCCTGGGCCAGATTCCAATG
CAAAGCAACCCATTAGAAAGACTCAACCACCTTGAAGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202856 representing NM_012279 Red=Cloning site Green=Tags(s) MEYPAPATVQAADGGAAGPYSSSELLEGQEPDGVRFDREARRLWEAVSGAQPVGREEVEHMIQKNQCLF TNTQCKVCCALLISESQKLAHYQSKKHANKVKRYLAIHGMETLKGETKKLSDQKSSRSKDKNQCCPICN MTFSSPVVAQSHYLGKTHAKNLKQSTKVEALHQNREMIIDPKFCSLCHATFNDPVMAQQHYVGKKHR KQETKLKLMARYGRLADPAVTDFFPAGKGYPCCKIVLNSIEQYQAHVSGFKHKNQSPKTVASSLGQIPM QRQPIQKDSTTLED TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_012279
ORF Size:	882 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_012279.2</u>
RefSeq Size:	3089 bp
RefSeq ORF:	885 bp
Locus ID:	23567
UniProt ID:	<u>Q9UL40</u>
Cytogenetics:	5q35.2
MW:	32.9 kDa