

Product datasheet for **RC202856L3V**

ZNF346 (NM_012279) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF346 (NM_012279) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF346
Synonyms:	JAZ; Zfp346
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012279
ORF Size:	882 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202856).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012279.2
RefSeq Size:	3089 bp
RefSeq ORF:	885 bp
Locus ID:	23567
UniProt ID:	Q9UL40
Cytogenetics:	5q35.2
Domains:	ZnF_U1, zf-C2H2
Protein Families:	Druggable Genome



[View online »](#)

MW: 32.9 kDa

Gene Summary: The protein encoded by this gene is a nucleolar, zinc finger protein that preferentially binds to double-stranded (ds) RNA or RNA/DNA hybrids, rather than DNA alone. Mutational studies indicate that the zinc finger domains are not only essential for dsRNA binding, but are also required for its nucleolar localization. The encoded protein may be involved in cell growth and survival. It plays a role in protecting neurons by inhibiting cell cycle re-entry via stimulation of p21 gene expression. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Apr 2015]