

## Product datasheet for **RC202885L3V**

### **LKB1 (STK11) (NM\_000455) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	LKB1 (STK11) (NM_000455) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LKB1
Synonyms:	hLKB1; LKB1; PJS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000455
ORF Size:	1299 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202885).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000455.4</a>
RefSeq Size:	3286 bp
RefSeq ORF:	1302 bp
Locus ID:	6794
UniProt ID:	<a href="#">Q15831</a>
Cytogenetics:	19p13.3
Domains:	pkinese, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase




[View online »](#)

**Protein Pathways:** Adipocytokine signaling pathway, mTOR signaling pathway

**MW:** 48.5 kDa

**Gene Summary:** This gene, which encodes a member of the serine/threonine kinase family, regulates cell polarity and functions as a tumor suppressor. Mutations in this gene have been associated with Peutz-Jeghers syndrome, an autosomal dominant disorder characterized by the growth of polyps in the gastrointestinal tract, pigmented macules on the skin and mouth, and other neoplasms. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized. [provided by RefSeq, Jul 2008]

### Product images:



[RC202885L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC202885L3V particle to overexpress human STK11-Myc-DDK fusion protein.