

Product datasheet for **RC202897L4V**

ATF3 (NM_001030287) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | ATF3 (NM_001030287) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | ATF3 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001030287 |
| ORF Size: | 543 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC202897). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001030287.2 |
| RefSeq Size: | 1935 bp |
| RefSeq ORF: | 546 bp |
| Locus ID: | 467 |
| UniProt ID: | P18847 |
| Cytogenetics: | 1q32.3 |
| Protein Families: | Transcription Factors |
| MW: | 20.6 kDa |



[View online »](#)

Gene Summary:

This gene encodes a member of the mammalian activation transcription factor/cAMP responsive element-binding (CREB) protein family of transcription factors. This gene is induced by a variety of signals, including many of those encountered by cancer cells, and is involved in the complex process of cellular stress response. Multiple transcript variants encoding different isoforms have been found for this gene. It is possible that alternative splicing of this gene may be physiologically important in the regulation of target genes. [provided by RefSeq, Apr 2011]