

Product datasheet for **RC202948A1V**

Human TACO1 (NM_016360) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TACO1 (NM_016360) AAV Particle
Tag:	Myc-DDK
Symbol:	TACO1
Synonyms:	CCDC44; MC4DN8
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202948 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGC**

ATGTCGGCTTGGGCTGCTGCCAGCCTAAGCAGGGCCGCTGCCCGATGCTTGGTGGCAGAGGCCCGGGG
TCAGGGCGGCTCCTCCGCGGACCCCGGCCCTCCACCCGAGCCCGGGGCTGCGGTGCCGCTCCGGG
CAGGACGCTGCACTTTACCGCGGCTGTCCCGCCGGGCACAACAAGTGGTCCAAAGTCAGGCACATCAAG
GGTCCGAAGGACGTCGAAAGGAGTCGCATCTTCTCAAACCTCTGTTGAACATCCGCCTGGCAGTGAAAG
AAGGAGGCCCAACCTGAGCACAACAGCAACCTGGCCAATATCTTAGAGGTGTGTCGCAGCAACATAT
GCCCAAGTCAACGATTGAGACGCACTGAAAATGGAGAAATCCAAGGACACTTATTTGCTGTATGAGGGT
CGAGGCCCTGGTGGCTCTTCTCTGCTCATCGAGGCATTATCTAACAGTAGCCACAAGTGCCAAGCAGACA
TTAGACATATCCTGAATAAGAATGGAGGAGTGATGGCTGTAGGAGCTCGTCACTCTTTTGACAAAAGGG
GGTGATTGTGGTTGAAGTGGAGGACAGAGAGAAGAAGGCTGTGAACCTAGAGCGTGCCCTGGAGATGGCA
ATCGAAGCAGGAGCTGAGGATGTCAAGGAACTGAAGATGAAGAAGAAAGAACGTTTTTAAATTTATTT
GTGATGCCTCTTCACTGCACCAAGTGAGGAAGAAGCTGGACTCCCTGGGCCTGTGTTCTGTCTCTGTGC
ACTAGAGTTTATCCCAACTCAAAGGTGCAGCTGGCTGAGCCCGACCTGGAACAGGCCGCACATCTCATT
CAGGCTCTCAGCAACCACGAGGATGTGATTCACGTCTATGATAACATTGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202948 protein sequence Red=Cloning site Green=Tags(s) MSAWAAASLSRAAARCLLARGPGVRAAPPRDPRPSHPEPRGCGAAPGRTLHFTAAVPAGHNKWSKVRHIK GPKDVERSRIFSCLKLNIRLAVKEGGPNPEHNSNLANILEVCRSKHMPKSTIETALKMEKSKDTYLLYEG RGPGGSSLLIEALSNSSHKCQADIRHILNKNGGVMAGARHSFDKKGVIIVVEVEDREKKAVNLERALEMA IEAGAEDVKETEDDEERNVFKFICDASSLHQVRKKLDSLGLCSVSCALEFIPNSKVQLAEPDLEQAAHLI QALSNHEDVIHVYDNIE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_016360
ORF Size:	891 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_016360.2</u>
RefSeq Size:	1479 bp
RefSeq ORF:	894 bp
Locus ID:	51204
UniProt ID:	<u>Q9BSH4</u>
Cytogenetics:	17q23.3
MW:	32.5 kDa