

## Product datasheet for **RC202951L3V**

### Syntaxin 6 (STX6) (NM\_005819) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Syntaxin 6 (STX6) (NM_005819) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Syntaxin 6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005819
ORF Size:	765 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202951).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005819.4</a>
RefSeq Size:	5074 bp
RefSeq ORF:	768 bp
Locus ID:	10228
UniProt ID:	<a href="#">O43752</a>
Cytogenetics:	1q25.3
Domains:	t_SNARE
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	SNARE interactions in vesicular transport



[View online »](#)

**MW:** 29.2 kDa

**Gene Summary:** Involved in intracellular vesicle trafficking.[UniProtKB/Swiss-Prot Function]