

Product datasheet for **RC202958L4V**

DNAJC30 (NM_032317) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNAJC30 (NM_032317) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNAJC30
Synonyms:	LHONAR; MC1DN38; WBSCR18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032317
ORF Size:	678 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202958).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032317.2
RefSeq Size:	2534 bp
RefSeq ORF:	681 bp
Locus ID:	84277
UniProt ID:	Q96LL9
Cytogenetics:	7q11.23
Domains:	DnaJ
Protein Families:	Transmembrane


[View online »](#)

MW: 26.1 kDa

Gene Summary: This intronless gene encodes a member of the DNAJ molecular chaperone homology domain-containing protein family. This gene is deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.23. [provided by RefSeq, Jul 2008]