

Product datasheet for **RC202985L3V**

Pygopus 2 (PYGO2) (NM_138300) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pygopus 2 (PYGO2) (NM_138300) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Pygopus 2
Synonyms:	1190004M21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_138300
ORF Size:	1218 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202985).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138300.3
RefSeq Size:	3227 bp
RefSeq ORF:	1221 bp
Locus ID:	90780
UniProt ID:	Q9BRQ0
Cytogenetics:	1q21.3
Protein Families:	Druggable Genome, Transcription Factors
MW:	41.1 kDa



[View online »](#)

Gene Summary: Involved in signal transduction through the Wnt pathway.[UniProtKB/Swiss-Prot Function]