

Product datasheet for **RC202994L4V**

PAAF1 (NM_025155) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PAAF1 (NM_025155) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PAAF1
Synonyms:	PAAF; Rpn14; WDR71
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025155
ORF Size:	1176 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202994).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025155.1
RefSeq Size:	1626 bp
RefSeq ORF:	1179 bp
Locus ID:	80227
UniProt ID:	Q9BRP4
Cytogenetics:	11q13.4
Domains:	WD40
MW:	42.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a WD repeat-containing protein involved in regulation of association of proteasome components. During HIV infection, the encoded protein is thought to promote provirus transcription through recruitment of the 19S regulatory complex. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jun 2012]