

## Product datasheet for **RC202996L3V**

### ETEA (FAF2) (NM\_014613) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ETEA (FAF2) (NM_014613) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ETEA
Synonyms:	ETEA; UBXD8; UBXN3B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014613
ORF Size:	1335 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202996).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014613.2</a>
RefSeq Size:	4515 bp
RefSeq ORF:	1338 bp
Locus ID:	23197
UniProt ID:	<a href="#">Q96CS3</a>
Cytogenetics:	5q35.2
Domains:	UBX, UAS
MW:	52.4 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is highly expressed in peripheral blood of patients with atopic dermatitis (AD), compared to normal individuals. It may play a role in regulating the resistance to apoptosis that is observed in T cells and eosinophils of AD patients. [provided by RefSeq, Jul 2008]