

## Product datasheet for **RC203014L1V**

### TM2D3 (NM\_078474) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TM2D3 (NM_078474) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TM2D3
Synonyms:	BLP2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_078474
ORF Size:	741 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203014).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_078474.2</a> , <a href="#">NP_510883.1</a>
RefSeq Size:	1421 bp
RefSeq ORF:	744 bp
Locus ID:	80213
UniProt ID:	<a href="#">Q9BRN9</a>
Cytogenetics:	15q26.3
Domains:	TM2
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

**MW:** 27.2 kDa

**Gene Summary:** The protein encoded by this gene contains a structural module related to that of the seven transmembrane domain G protein-coupled receptor superfamily. This protein has sequence and structural similarities to the beta-amyloid binding protein (BBP), but, unlike BBP, it does not regulate a response to beta-amyloid peptide. This protein may have regulatory roles in cell death or proliferation signal cascades. Several alternatively spliced transcript variants of this gene are described but the full length nature of some variants has not been determined. Multiple polyadenylation sites have been found in this gene. [provided by RefSeq, Jul 2008]