

Product datasheet for **RC203044L3V**

NMD3 (NM_015938) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NMD3 (NM_015938) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NMD3
Synonyms:	CGI-07
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015938
ORF Size:	1509 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203044).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015938.2
RefSeq Size:	2758 bp
RefSeq ORF:	1512 bp
Locus ID:	51068
UniProt ID:	Q96D46
Cytogenetics:	3q26.1
Domains:	NMD3
MW:	57.6 kDa


[View online »](#)

Gene Summary:

Ribosomal 40S and 60S subunits associate in the nucleolus and are exported to the cytoplasm. The protein encoded by this gene is involved in the passage of the 60S subunit through the nuclear pore complex and into the cytoplasm. Several transcript variants exist for this gene, but the full-length nature of only two have been described to date. [provided by RefSeq, Feb 2016]