

Product datasheet for **RC203055L4V**

MICS1 (GHITM) (NM_014394) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MICS1 (GHITM) (NM_014394) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MICS1
Synonyms:	DERP2; HSPC282; MICS1; My021; PTD010; TMBIM5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014394
ORF Size:	1035 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203055).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014394.2
RefSeq Size:	2449 bp
RefSeq ORF:	1038 bp
Locus ID:	27069
UniProt ID:	Q9H3K2
Cytogenetics:	10q23.1
Protein Families:	Transmembrane
MW:	37.2 kDa



[View online »](#)

Gene Summary:

Required for the mitochondrial tubular network and cristae organization. Involved in apoptotic release of cytochrome c.[UniProtKB/Swiss-Prot Function]