

Product datasheet for **RC203061A1V**

Human BCCIP (NM_078468) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human BCCIP (NM_078468) AAV Particle
Tag:	Myc-DDK
Symbol:	BCCIP
Synonyms:	TOK-1; TOK1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203061 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGTCCAGGTCTAAGCGGCGTGCCGTGGAAGTGGGGTCCGCAGCCGCCGGATCCCCAGTCCAGC
GCGACGAGGAAGAGGAAAAAGAAGTCGAAAATGAGGATGAAGACGATGATGACAGTGACAAGGAAAAAGGA
TGAAGAGGACGAGGTCATTGACGAGGAAGTGAATATTGAATTTGAAGCTTATTCCTATCAGATAATGAT
TATGACGGAATTAAGAAATTACTGCAGCAGCTTTTCTAAAGGCTCCTGTGAACACTGCAGAATAACAG
ATCTCTTAATTCAACAGAACCATTGGGAGTGTGATTAAGCAAACGGATGTTTCAGAAGACAGCAATGA
TGATATGGATGAAGATGAGGTTTTGGTTTCATAAGCCTTTTAAATTTAACTGAAAGAAAGGGTACCCAG
TGTGTTGAACAAATTCAAGAGTTGGTTCTACGCTTCTGTGAGAAGAACTGTGAAAAGAGCATGGTTGAAC
AGCTGGACAAGTTTTTAAATGACACCACCAAGCCTGTGGGCCTTCTCCTAAGTAAAAGATTCTTAATGT
CCCTCCACAGATCGCTCTGCCCATGTACCAGCAGCTTCAGAAAGAACTGGCGGGGACACAGAACCAAT
AAGCCATGTGGGAAGTGCTACTTTTACCTTCTGATTAGTAAGACATTTGTGAAGCAGAAAAACAATT
CCAAAAAGAACTAGCAACAAAAAGAAAGCTGCGTTAATGTTTGCAATGCAGAGGAAGATTTTCTA
TGAGAAGGCAATTCTCAAGTTCAACTACTCAGTGCAGGAGGAGAGCGACACTTGTCTGGGAGGCAATGG
TCTTTTGATGACGTACCAATGACGCCCTTGCGAACTGTGATGTTAATTCCAGGCGACAAGATGAACGAAA
TCATGGATAAACTGAAAGAATATCTATCTGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203061 protein sequence Red=Cloning site Green=Tags(s) MASRSKRRAVESGVPQPPDPVQRDEEEKEVENEDDDSDKEKDEEDEVIDEEVNIEFEAYSLSDND YDGIKKLLQQLFLKAPVNTAELTDLLIQNHIGSVIKQTDVSEDSNDDMDEDEVFGFISLLNLTERKGTQ CVEQIQELVLRFCCKNCEKSMVEQLDKFLNDTTKPVGLLLSERFINVPPQIALPMYQQLQKELAGAHRTN KPCGKCYFYLLISKTFVEAEKNNSKKKPSNKKKAALMFANAEFFFEYKAILKFNYSVQEEESDTCLGGKW SFDDVPMTPLRTVMLIPGDKMNEIMDKLKEYLSV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_078468
ORF Size:	942 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_078468.1</u>
RefSeq Size:	1261 bp
RefSeq ORF:	945 bp
Locus ID:	56647
UniProt ID:	<u>Q9P287</u>
Cytogenetics:	10q26.2
MW:	36.1 kDa